

BRAND	MODEL	PLATES			
		PCRL-096-100	PCRL-Q96-100	PCRL-L96-100	PCRL-W96-100
Applied Biosystems (Life Technology)	7500Fast			•	
	7900HT Fast Model (0,1 ml)			•	
	StepOne Plus			•	
	ViiA7 Fast Model (0,1ml)			•	
	QuantStudio Fast Model (0,1ml)			•	
	QuantStudio 3/5 Fast (0,1 ml)			•	
	5700/7700		•		
	3100/3130 (xL) (0,1 ml Fast 96 well)			•	
	3100/3130 (xL) (0,2 ml 96 well)		•		
	3700/3730 (xL) (0,1 ml Fast 96 well)			•	
	3700/3730 (xL) (0,2 ml 96 well)		•		
	3500 (xL) (0,1 ml Fast 96 well)			•	
3500 (xL) (0,2 ml 96 well)		•			
Roche	LightCycler 480 & 480II				•
	LightCycler 96				•
Bio-Rad	iCycler iQ	•	•		
	iQ5	•	•		
	MyiQ	•	•		
	Chromo4	•	•		
	Opticon I & II				•
	CFX96				•
Stratagene/IT-IS /Eppendorf/ Thermo/Biometra	Mx4000/3000P/3005P	•	•		
	AriaMx				•
	MasterCycle ep Realplex	•	•		
	Toptical thermocycler	•	•		